Peregrine

Manufactured by:

Peregrine Surgical Ltd.

51 Britain Drive

New Britain, PA 18901

Tel: 215-348-0456

Fax: 215-348-5526

E-mail: oph@peregrine-surgical.com

Website: www.peregrinesurgical.com

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EN

Indications for Use

For photocoagulation and illumination during ophthalmic surgery. This device delivers illumination as well as laser energy to target tissue, causing coagulation. Spot size can be varied by altering the distance between the tissue and the probe tip.

Instructions for Use

The Laser Probe can only be used in combination with a medical laser in the Vis-NIR transmission range of 532nm to 810nm with a maximum power output of 2.5 Watts and that is appropriate for photocoagulation. The light portion of the device is designed for use with the Alcon Constellation Vision System. The Constellation uses a 50 Watt Xenon arc lamp that emits light with an output of 30 lumens within the wavelength range of 400nm to 700nm. Please verify that the Laser Probe is compatible with the ophthalmic system with which it is used. See laser operations manual for safety precautions and to ensure both the laser and light source meet above requirements; if the requirements are not met do not use this device. The device is intended to transmit energy; but does not generate energy. Optical radiation testing resides with the manufacturer of the laser.

Packaging

The Laser Probes are supplied sterile. The Laser Probes are ethylene oxide sterilized and guaranteed to be sterile until the expiration date stated on the product label. Each Laser Probe is packaged in a double peel pouch for ease of introduction into the sterile field.

Preparation

- Please visually inspect the packaging of the sterile Laser Probe for possible damages and expiration date. Probes from previously opened or damaged packaging must be deemed unsterile.
- Unpack the Laser Probe in a sterile environment.
- Set the medical laser equipment in a safe setting following the instructions of the system manual.
- Remove Laser Probe from protections tray, remove band and uncoil.
- The proximal end (SMA connector) is inserted into the receptor port of the medical laser equipment following the instructions of the system manual.
- Insert the illumination connector into the appropriate light source or light source adapter.
- The distal end handpiece is used as indicated.

Performance testing

To test the performance for the safe use of the Laser Probe, the image of the aiming beam can be inspected before the start of the laser procedure. For this purpose the tip of the probe is pointed straight to a bright (white) surface at a distance of about 20cm (8"). The usual colored aiming beam is activated on the medical laser equipment.

Laser Probe

When the image is circular and evenly illuminated, the performance of the Laser Probe is acceptable. If the image is faint or unequal, the Laser Probe should not be used.

Laser Probe

Each probe has been tested to ensure efficient transmission of light energy at the operational wavelengths. The average transmission efficiency for each gauge is as follows: 20ga – 96%, 23ga – 60%, 25ga – 38%.

Directions for Use

Insert probe through incision. Positioning the Laser Probe closer to its intended target will provide a decreased laser spot size as well as a focused light field. Backing probe away from its intended target will result in a larger laser spot size and wider illuminating field.

Laser Probe

Warning
Avoid unnecessary illumination of the retina using the aiming beam. Please read carefully the instructions for full clinical use provided by the manufacturer of the medical laser equipment. Do not apply excessive stress to the Laser Probe as to avoid product damage. Avoid contact of the Laser Probe tip with human tissue, blood or other instruments, because of the risk of uncontrolled scatter or damage to the Laser Probe. Please store the product in a dry and clean environment. Reuse may not guarantee optimal performance, safety or sterility and may cause patient injury or infection.

Laser Probe

The use of flammable anesthetics or oxidizing gases such as nitrous oxide and oxygen should be avoided. Some materials, for example cotton or wool when saturated with oxygen, may be ignited by the high temperatures produced in normal use of the laser equipment. The solvents of adhesives and flammable solutions used for cleaning and disinfecting should be allowed to evaporate before the laser equipment is used. Attention should also be drawn to the danger of ignition of endogenous gases.

Laser Probe

User to report any serious incident that has occurred in relation to this device to the manufacturer and the authority having jurisdiction in the incident locale.

CAUTION

(USA) Federal law restricts this device to sale by or on the order of a licensed physician.

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