 **Manufactured by:**

**Peregrine Surgical Ltd.**  
**51 Britain Drive**  
**New Britain, PA 18901**  
**Tel: 215-348-0456**  
**Fax: 215-348-5526**  
**E-mail: [oph@peregrine-surgical.com](mailto:oph@peregrine-surgical.com)**  
**Website: [www.peregrinesurgical.com](http://www.peregrinesurgical.com)**

**GB Care & Handling instructions**  
**F Instructions d'entretien et de manipulation**  
**D Gebrauchs- und Pflegeanweisungen**  
**E Instrucciones para su cuidado y manipulación**  
**P Instruções de Cuidado & Manuseamento**  
**I Istruzioni per la cura e la manipolazione**  
**S Skötselinstruktioner**

 : Medical Device

 : Prescription Only

 : Sterilized using ethylene oxide

 : Caution

 : Do not re-use

 : Do not use if package is damaged

 : Do not resterilize

 : Does not contain latex

 : Consult instructions for use

 : Use by Date

 : Batch Code

 : Keep dry

 : Keep away from sunlight

 : Double Sterile Barrier System

PK401N02 Rev. 1

 :

ADVENA LTD.

Tower Business Centre, 2<sup>nd</sup> Fl.

Tower Street, Swatar, BKR 4013

Malta

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## EN

**Indications for Use Statement:** For photocoagulation and aspiration during ophthalmic surgery.

### Instructions for Use

The Laser Probe can only be used in combination with a medical laser in the Vis-NIR transmission range of 532 nm to 810 nm. Please verify that the Laser Probe is compatible with the ophthalmic system with which it is used. The device is intended to transmit energy but does not generate energy. Optical radiation testing resides with the manufacturer of the laser. See laser operations manual for safety precautions and to assure that the laser meets the above requirements.

### Packaging

The Laser Probes are supplied sterile. The Laser Probes are ethylene oxide sterilized and guaranteed to be sterile until the expiration date stated on the product label. Each Laser Probe is packaged in a double peel pouch for ease of introduction into the sterile field.

### Preparation

- Please visually inspect the packaging of the sterile Laser Probe for possible damages and expiration date. Probes from previously opened or damaged packaging must be deemed unsterile.
- Unpack the Laser Probe in a sterile environment.
- Set the medical laser equipment in a safe setting following the instructions of the system manual.
- Remove Laser Probe from protections tray, remove band and uncoil.
- The proximal end (SMA connector) is inserted into the receptor port of the medical laser equipment following the instructions of the system manual.
- Attach luer to the aspiration line connected to the operating system.
- The distal end handpiece is used as indicated.

### Performance testing

To test the performance for the safe use of the Laser Probe, the image of the aiming beam can be inspected before the start of the laser procedure. For this purpose the tip of the probe is pointed straight to a bright (white) surface at a distance of about 20cm (8"). The usual colored aiming beam is activated on the medical laser equipment.

When the image is circular and consistent, the performance of the Laser Probe is acceptable. If the image is faint or unequal, the Laser Probe should not be used.

### Directions for Use

Insert probe through incision. Positioning the Laser Probe closer to its intended target will provide a decreased laser spot size. Backing probe away from its intended target will result in a larger laser spot size.

Activate aspiration, through control on operating system, when necessary. In case of tip occlusion, swiftly compress reflux button in order to free tip of unwanted material.

### Warning

Avoid unnecessary exposure of the retina using the aiming beam.

Please read carefully the instructions for full clinical use provided by the manufacturer of the medical laser equipment.

Do not apply excessive stress to the Laser Probe as to avoid product damage.

Avoid contact of the Laser Probe tip with other instruments, because of the risk of uncontrolled scatter or damage to the Laser Probe.

Please store the product in a dry and clean environment.

Reuse may not guarantee optimal performance, safety or sterility and may cause patient injury or infection.

The use of flammable anesthetics or oxidizing gases such as nitrous oxide and oxygen should be avoided. Some materials, for example cotton or wool when saturated with oxygen, may be ignited by the high temperatures produced in normal use of the laser equipment. The solvents of adhesives and flammable solutions used for cleaning and disinfecting should be allowed to evaporate before the laser equipment is used. Attention should also be drawn to the danger of ignition of endogenous gases.

User to report any serious incident that has occurred in relation to this device to the manufacturer and the authority having jurisdiction in the incident locale.

### Caution

(USA) Federal law restricts this device to sale by or on the order of a licensed physician.

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